



### Challenge

A simple, fast, and sensitive detection of total fluorine in wastewater

### Solution

With regard to the determination of a sum parameter of fluorinated compounds, a species-unspecific response of high-resolution continuum source graphite furnace molecular absorption spectrometry (HR-CS GF MAS) was developed.

## A Novel Method for the Fast, Sensitive, and Simple Analysis of Total Fluorine in Wastewater

### Introduction

The introduction of carbon-fluorine bonds into organic compounds is known to profoundly influence their chemical and physical properties when compared to their non-fluorinated analogues<sup>1</sup>. Many of the most important products of the chemical and life-science industries, from pharmaceuticals, through fine and specialty chemicals, to polymers, rely on organic fluorine chemistry for their useful properties. The accelerating production volumes of fluorinated organic substances expectably lead to their increasing release into the environment. Some of these substances are already listed as persistent organic pollutants by the Stockholm Convention (Annex B -restriction) and the European water framework directive (WFD).

Per- and polyfluoroalkyl substances (PFAS) are a large class of synthetic chemicals that present numerous analytical challenges, including their widespread presence in a variety of environmental samples<sup>2</sup>. In the US, PFAS may contaminate public drinking-water systems that

serve an estimated 19 million people<sup>3</sup>. The United States Environmental Protection Agency (EPA) is continuing to aggressively implement their PFAS Action Plan – the most comprehensive cross-agency plan ever to address an emerging chemical of concern<sup>2</sup>. The EPA has established a non-enforceable health advisory level of 70 parts per trillion (ppt) for the sum of pentadecafluorooctanoic acid (PFOA) and heptadecafluorooctanesulfonic acid (PFOS) for drinking water. The EPA is developing a potential rapid screening tool to identify total PFAS presence and absence in 2021. This eventual standard operating procedure will be used to quantify Total Organic Fluorine (TOF)<sup>4</sup>.

Therefore, a sensitive, fast, and simple TOF detection method is needed for environment pollution monitoring and control. The method is based on a combination of solid-phase extraction (SPE) with high-resolution-continuum source graphite furnace molecular absorption spectrometry (HR-CS

GF MAS). In another application note, "Determination of Extractable Organically Bound Fluorine (EOF) in Surface Water", the solid phase extraction (SPE) procedure has been fully discussed for the extraction of fluorinated compounds. HR-CS GF MAS was used to detect total fluorine via in-situ

formation of gallium(I) fluoride (GaF) molecules. In this study, we focus on improving the performance of Total Fluorine (TF) analysis in wastewater by optimizing HR-CS GF MAS method and applying novel calibration strategies.

## Method Development Strategy

For the quantification of fluorinated organic compounds as a sum parameter, a species-unspecific response for HR-CS GF MAS method needs to be developed and optimized. From our past research, gallium has proven to enable the most sensitive fluorine detection via diatomic GaF formation. The characteristic molecular absorption of the GaF is then detected by HR-CS-GF MAS. To achieve the best signal, the furnace tube was conditioned with the molecule forming agent (Ga) and modifiers (Pd/Mg/Zr/Ba) before the sample was injected. From the literature, it's known that

melting point/vapor pressure of the fluorinated compounds significantly impacts the recovery rates of the different compounds. Thus, the drying and pyrolysis procedure in the furnace program were optimized to keep possible losses of fluorinated compounds as low as possible. To evaluate the optimization approach, 24 different fluorine compounds were mixed for the QC sample. These 24 compounds were selected with different numbers of fluorine atoms and boiling points. Calibration strategies were developed to obtain the best recovery rate for this QC sample.

## Materials and Methods

### Reagents

- Ultrapure grade HNO<sub>3</sub>
- 0.05% TritonX-100
- Pd/Mg/Zr modifier (1 g/L Pd, 0.5 g/L Mg, 0.02 g/L Zr)
- Ba modifier (10 mg/L)
- Ga solution (10 g/L)
- Zr stock solution (1 g/L)
- Certificated F stock solution (0.1 g/L F ICP standard as NaF)
- Heptadecafluorooctanesulfonic acid solution (PFOS, 100 µg/mL in methanol, analytical standard)
- Pentadecafluorooctanoic acid solution (PFOA, 100 µg/mL in methanol, analytical standard)

### Samples

- C sample: 24 fluorinated organic compounds mixture (TOF about 128 µg/L)
- Wastewater sample

### Instrumentation

The contrAA 800 G graphite furnace AAS with AS-GF autosampler controlled by ASpect CS software was used for

the determination of fluorine. Before use, the graphite tube was first coated with zirconium (35 µL stock solution, 6 times). The fluorine content was determined by measuring the intensity of the gallium monofluoride absorption band. During the analysis, fluorine is converted stoichiometrically to GaF in the presence of gallium. To improve the sensitivity, a conditioning procedure was included in the furnace program. Each sample was analyzed with three replicates. Table 2 and Table 3 show the method settings and furnace program.

Table 1: Instrument specifications

Parameter	Specification
Instrument	contrAA 800 G
Furnace tube type	PIN platform
Autosampler	AS-GF
Injected volume	4–20 µL (standards), 20 µL (sample)
Rinsing solution	2% HNO <sub>3</sub> , 0.05% TritonX-100

Table 2: Method settings and evaluation parameters

Molecule	Wave-length [nm]	No. of eval. pixels	T Pyr. [°C]	T Atomis. [°C]	Ramp [°C/s]	Meas. Time [s]	Modifier	Baseline correction
GaF	211.248	5	500	1,450	1,200	5	5 µL Pd/Mg/Zr modifier* 5 µL Ga solution* 3 µL Ga solution 5 µL Ba solution	IBC

Note: \*-use in the conditioning step

Table 3: Furnace program for detection of GaF molecule

Step	Name	Temp (°C)	Ramp (°C/s)	Hold (s)	Gas purge
1	Drying*	80	6	5	max
2	Drying	100	6	5	max
3	Drying	160	10	5	max
4	Drying	350	25	10	max
5	Drying	1100	500	12	max
6	Drying	70	NP	10	max
7	Drying^	80	2	10	max
8	Drying	110	5	15	max
9	Pyrolysis	250	50	5	max
10	Pyrolysis	500	200	5	max
11	Gas adaption	500	0	5	stop
12	Atomization	1,450	1,200	5	stop
13	Clean	2,450	1,200	4	max

Note: \*-modifier injection for pretreatment, ^-sample injection

### Calibration Strategy

From the past research and literature, the melting point/vapor pressure of the fluorinated compounds plays an important role for total fluorine analysis. If the vapor pressure at a given temperature of the respective substance is high, the analyte is transferred into the gas phase and carried away by the argon stream during the furnace program before detection. Thus, the calibration standards should be selected to have similar behavior to the samples. To find the best calibration strategy for the waste water, three different calibration sets were prepared. Table 4 shows the components of each calibration set.

For the inorganic fluoride calibration (inorganic F), a 1,000 ppm NaF stock solution was used for preparing calibration standards. Three common fluorinated compounds, PFOS, PFOA and HFPO-DA, were selected to prepare organic

fluorine (organic F) standards. The inorganic and organic mixture standard was also prepared with NaF, PFOS, PFOA and HFPO-DA. Table 5 shows the concentrations and absorbances of the standards, as well as the linearity of calibration. Based on the calibration curve and measured blank values, the lowest LOD of 4.00 µg/L was achieved with the inorganic and organic mixture calibration. In Table 5, the inorganic fluoride calibration standards show the highest signal while the organic calibration standards show the lowest. The different signal responses of organic and inorganic substances depends on their thermal stability. Volatile PFAS are partially lost during the drying and pyrolysis steps of the furnace program. However, inorganic fluoride is thermally stable during these steps. In Table 6, the clear and interference-free GaF absorbance spectra (blue) at 211.248 nm is shown from the three calibration sets.


Table 4: Concentration sets

Calibration Set	F Compounds	Ratio
Inorganic F	NaF	
Organic F	PFOS, PFOA, HFPO-DA	1:1:1
Inorganic and organic F mixture	NaF, PFOS, PFOA, HFPO-DA	1:1:1:1

Table 5: Type and quality of calibration

Name	Conc.	Inorganic F	Organic F	Inorganic and organic mixture
Unit	µg/L	Abs	Abs	Abs
Cal-Zero1	0	0.1726	0.1153	0.1709
Cal-Std1	125	1.6650	0.5177	1.1601
Cal-Std2	250	3.0402	1.0367	1.7833
Cal-Std3	375	4.4542	1.4847	2.4628
Cal-Std4	500	5.9955	1.9403	3.0694
R <sup>2</sup>		0.9993	0.9999	0.9808
LOD (µg/L)		7.74	12.58	4.00

Table 6: Characteristic signal shape and spectral vicinity of the analyte line



\*blue: GaF analyte signal, red: background signal

## Results and Discussion

### 1. Result of inorganic F calibration strategy

With the inorganic F calibration strategy, the QC standard of 24 fluorinated organic compounds shows a 38% recovery rate (Table 7). Differences in signal response result in the low recovery rate for PFAS when using NaF solutions as the calibration standard. When the QC sample is spiked with 100 µg/L F (as NaF), it shows a good spike recovery rate of 114%.

Table 7: Measurement results by inorganic F calibration

Samples	Conc. (µg/L)	Recovery (%)
Blank water	<LOQ	
QC sample	48.9	38
QC sample spiked with 100 µg/L (NaF)	162.5	114
Wastewater	348.2	

### 2. Result of organic F calibration strategy

With the organic F calibration strategy, the QC sample's recovery is 190%. This overestimation of PFAS is due to the different thermal stabilities of the various organic fluorinated compounds. When the QC is spiked with 100 µg/L NaF, a 275% spike recovery rate is obtained due to the stronger stability and signal response of NaF. Both the QC and wastewater samples were spiked with 125 µg/L F from the TOF mixture. This TOF mixture is composed of the same compounds as the organic calibration standards. The spike recovery rates for both the QC and wastewater are 74%.

Table 8: Measurement results by organic F calibration

Samples	Conc. (µg/L)	Recovery (%)
Blank water	<LOQ	
QC sample	224.7	174
QC sample spiked with 125 µg/L F (TOF*)	316.7	74
QC sample spiked with 100 µg/L F (NaF)	628.1	275
Wastewater	4882.5	
Wastewater (1:50 dilution) spiked with 125 µg/L F (TOF*)		74

\*TOF = PFOS, PFOA, HFPO-DA (ratio 1:1:1)

### 3. Result of inorganic and organic mixture calibration strategy

From Table 9, the organic and inorganic F mixture calibration strategy provides the best approximation of sample analyte characteristics. The QC samples recovery rate is 100%. The QC sample and wastewater sample were spiked with 125 µg/L F from the TF mixture. The compounds in the TF mixture are the same as the inorganic and organic mixture calibration standards. The spike recovery rates are 85% and 86% for the QC sample and wastewater sample, respectively. The QC is also spiked with inorganic F and the recovery rate is overestimated as 167%. This is also due to the higher thermal stability of inorganic NaF. Comparison between the three calibration strategies demonstrates the dramatic influence the calibration components can impact the analysis results. For this specific QC sample, the inorganic and organic F mixture calibration strategy is the optimal choice.

Table 9: Measurement results by inorganic and organic mixture calibration

Samples	Conc. ( $\mu\text{g/L}$ )	Recovery (%)
Blank water	<LOQ	
QC sample	127.6	100
QC sample spike with 125 $\mu\text{g/L}$ F (TF*)	252.6	85
QC standard spike with 100 $\mu\text{g/L}$ F (NaF)	380.9	167
Wastewater	1936.5	
Wastewater (1:50) spike with 125 $\mu\text{g/L}$ F (TF*)		86

\*TF = NaF, PFOS, PFOA, HFPO-DA (ratio 1:1:1:1)

#### 4. Long-term stability test

To evaluate the long-term stability of this method, the QC sample was tested on three separate days each one week apart. The three tests all used Test 1's calibration curve without recalibration. In this period, the RSD is 5.7% of those three tests (Table 10). This demonstrates the excellent long-term stability of the method.

Table 10: Long-term stability test

	Conc. ( $\mu\text{g/L}$ )	RSD (5 Replicates)
Test 1 (week 1)	127.6	3.8 %
Test 2 (week 2)	113.8	3.0%
Test 3 (week 3)	120.2	2.7%
Average Conc.	120.5	
RSD	5.7%	

#### 5. Blank value

During method optimization, relatively high blank values were observed. The carry-over effect in the sample introduction system was considered to be the largest contributing factor. By adding  $\text{HNO}_3$  and TritonX-100 in the rinsing solution, the blank value was lower. It's very important to use F-free modifiers. Using high purity salt, such as gallium (III) nitrate hydrate, to prepare modifiers is recommended.

## Summary

In this study, a fast, easy, and sensitive complementary method for total fluorine analysis in wastewater is presented. This study shows that the MAS method can be successfully applied using the contraAA 800 G for the measurements. The optimized furnace program and calibration strategy provide high sample throughput, sensitivity, and accuracy. Only 3.5 minutes/replicate is needed for each sample. 100% recovery rate of QC sample is successfully achieved with an inorganic and organic mixture calibration strategy. The lowest LOD of 4 ppb is also achieved the inorganic and organic mixture calibration strategy. The results of fluorine determination by HR-CS-GF MAS show good reproducibility and long-term stability. Sample dilution and spiking can be easily achieved with an AS-GF autosampler. No additional cleaning steps are required for this method.



Figure 1: contraAA 800 G

## References

- [1] Richard, D.; Chambers, F.R.S.; FLUORINE IN ORGANIC CHEMISTRY. Blackwell Publishing Ltd. 2004
- [2] EPA Press Office; AGGRESSIVELY ADDRESSING PFAS AT EPA; <https://www.epa.gov/newsreleases/aggressively-addressing-pfas-epa>; 2020 Jan
- [3] Cheryl H.; WHY LIMITING PFAS IN DRINKING WATER IS A CHALLENGE IN THE US. Chemicals and Engineer News. 2020 Jul; 98 (27)
- [4] <https://www.epa.gov/water-research/pfas-analytical-methods-development-and-sampling-research>
- [5] Metzger, M.; Ley, P.; Sturm, M. and Meermann, B.; SCREENING METHOD FOR EXTRACTABLE ORGANICALLY BOUND FLUORINE (EOF) IN RIVER WATER SAMPLES BY MEANS OF HIGH-RESOLUTION-CONTINUUM SOURCE GRAPHITE FURNACE MOLECULAR ABSORPTION SPECTROMETRY (HR-CS GF MAS). Anal Bioanal Chem. 2019 Jul; 411 (19): 4647-4660

This document is true and correct at the time of publication; the information within is subject to change. Other documents may supersede this document, including technical modifications and corrections.

Trademark notice: The brand names of the third-party products specified in the application note are usually registered trademarks of the respective companies or organizations.

## Headquarters

Analytik Jena GmbH+Co. KG    Phone +49 3641 77 70    info@analytik-jena.com  
Konrad-Zuse-Strasse 1    Fax +49 3641 77 9279    www.analytik-jena.com  
07745 Jena · Germany

Version 1.1- Author: SS  
en · 04/2025  
© Analytik Jena GmbH+Co. KG | Pictures: pexels